



ClearSounds®

Hear how good life can be

CS-CLA7BT

Bluetooth Amplified Neckloop Accessory



ENGLISH

ESPAÑOL

FRANÇAIS

Congratulations on your purchase of the CLA7BT Amplified Powered Neckloop from ClearSounds® Communications, Inc. The CLA7BT has been developed to facilitate the use of a Bluetooth capable cellular phone, landline phone, or audio device with a t-coil equipped hearing aid for hands-free wireless operation.

Features:

- 30 dB Adjustable Amplification
- Compatible with standard Bluetooth Version 1.2/2.1 Devices
- 30 foot range (maximum)
- User-friendly controls and easy pairing
- Great for mobile, office and home applications.

The CLA7BT generates an amplified induction (loop) signal that is picked up by the telecoil of your hearing aid(s). This option in the HA is activated by setting the hearing aid(s) to the "T" or "MT" position. Please read this User Manual carefully in order to get the very best results from your CLA7BT. This will tell you how to set up the CLA7BT correctly and explain all the features to you.

The CLA7BT is compatible with most cellular phones and audio devices which support Bluetooth wireless technology. We hope it will

give you many years of successful operation. Please keep this manual for easy reference.

Please visit our web site at www.clearsounds.com to:

- Register your product
- Download additional manuals and troubleshooting tips

PLEASE SAVE THESE INSTRUCTIONS

IMPORTANT SAFETY INSTRUCTIONS

This product is a mobile solution built on Bluetooth technology. Please read the instructions before using your CLA7BT.

Electromagnetic Interference/Compatibility: Nearly every electronic device is susceptible to electromagnetic interference (EMI) if inadequately shielded, designed, or otherwise configured for electromagnetic compatibility.

Facilities: To avoid electromagnetic interference and/or compatibility conflicts, turn off your CLA7BT in any facility where posted notices instruct you to do so. Hospitals or health care facilities may be using equipment that is sensitive to external RF energy.

Aircraft: When instructed to do so, turn off your device when on board an aircraft. Any use of a

device must be in accordance with applicable regulations per airline crew instructions.

Medical Devices: If you use any personal medical devices, consult the manufacturer of your device to determine if it is adequately shielded from RF energy. Your physician may be able to assist you in obtaining this information.

WARNING: *To avoid the possibility of electrical shock, do not use this product near water (For example, near a bathtub, wash basin, kitchen sink or laundry tub, in a wet basement or near a swimming pool) while in the bathtub or shower or when you are wet.*

WARNING: *Avoid using the telephone during electrical storms. There may be a remote risk of electric shock from lightning. If the electrical storm is present in your immediate area, possibility of electrical shock exists.*

WARNING: *Do not use a phone if you suspect a gas leak in the area. The phone's electrical contacts could create a spark, which could ignite any heavy concentration of leaking gas.*

Table of Contents

Safety Precautions and Instructions.....	3
Feature Location	5
Important Facts on Charging the CLA7BT	6
Pairing your CLA7BT to your Cell Phone.....	8

Pairing your CLA7BT to Additional Devices ...	10
Wearing the Unit.....	10
Turning your CLA7BT On and Off	11
Making and Receiving Calls Using Bluetooth	12
LED Indicator	16
Audio Tones	16
Using your CLA7BT as a Personal Audio	
Neckloop.....	17
Troubleshooting Guide.....	18
Warranty Information.....	20

FEATURE LOCATION



