



QH2 BLUETOOTH
PHONE HUB SYSTEM
USER'S GUIDE

Welcome!

Thank you for purchasing the ClearSounds QH2 Bluetooth Phone Hub System.

We hope that you enjoy the hands-free conversation that your new system provides. The hub amplifies and sends sound between phones and Bluetooth headsets, streamers or neckloops.

Some of the great features you'll enjoy include:

- **Hands-free conversation within 30 feet of the hub**
- **Adjustable tone and volume settings**
- **Compatible with Bluetooth listening devices (headsets, streamers or neckloops)**
- **Compatible with most telephones with keys on the base and with all single line cordless phones**

Contact information

Please contact us with any questions that you might have. We are happy to assist you!

ClearSounds Communications, Inc.
1743 Quincy Avenue, Suite 143 | Naperville, IL 60540 USA

800-965-9043 (toll-free) | www.clearsounds.com

Safety Precautions

Carefully read and observe the warnings and cautions in this manual and on the equipment.



Warning: Warnings must be observed to prevent bodily injury.

WARNING: Use the equipment only in the described manner to avoid bodily injury or damage to the equipment.

WARNING: Locate and use the hub away from water and damp areas, such as swimming pools, bathtubs, sinks and damp basements, and do not use the hub when you are wet, to avoid electrical shock. If the hub is submerged, unplug the phone base from AC power before you retrieve it.

WARNING: Do not overload extension cords or power strips. This can result in electric shock.

WARNING: Use this hub only with the supplied AC power adapter assembly. Other types of power adapter assemblies can damage the hub and cause personal injury.

WARNING: Avoid using the hub with a corded phone during an electrical storm. There is a slight risk of electrical shock from lightning during a storm.

WARNING: The hub can amplify volume that can cause hearing loss. For each user, adjust the volume to personal preferences.

WARNING: Do not spill liquid on the hub or clean it while it is plugged in to AC power. This can cause a short circuit, a fire or electric shock.

WARNING: Do not push any type of object into the hub. Contact with dangerous voltage points can occur, resulting in a fire.

Safety Precautions



Caution: Cautions indicate practices that could harm the hub or other objects.

CAUTION: Do not place objects on the cables or cords, and place them away from foot traffic.

CAUTION: Locate the hub away from heat sources of any type. Keep slots and openings on the components uncovered to enable ventilation and protection against overheating.

CAUTION: Connect the supplied AC power adapter assembly into a surge protector. Power surges can damage the equipment.

CAUTION: Use only a damp cloth to clean the unit. Do not use liquid or aerosol cleaners.

CAUTION: In these circumstances, unplug the hub from AC power, and refer to qualified service personnel:

- Liquid spills in the hub, or it is exposed to rain or precipitation.
- The hub does not operate according to descriptions in this manual, or performance distinctly changes
- The hub is dropped or the casing is damaged.

Contents

Contents	
Setup	5
Unpack	5
Corded Phone Setup.....	6
Connect the Hub as a Handset Amplifier	7
Check for Compatibility.....	8
Bluetooth Pairing and Disconnection.....	9
Answer or Place a Call Using a Corded Phone System	11
Cordless Phone Setup	14
Connect the Hub.....	15
Bluetooth Pairing and Disconnection.....	16
Answer or Place a Call Using a Cordless Phone System	17
Troubleshooting	20
Compliance.....	21
Specifications	24
Warranty.....	25
Accessories	26

Setup

Unpack

Ensure that your packing box includes all of these parts. If any parts are missing or damaged, please contact your vendor or ClearSounds.

IMPORTANT: YOUR RECEIPT IS REQUIRED FOR WARRANTY. PLEASE RETAIN YOUR RECEIPT.



Hub



Phone line splitter and
Hub-to-splitter cable




AC adapter

For Use with Corded Phones



Hub-to-base cable

Corded Phone Setup



Corded Phone Setup

Connect the Hub as a Handset Amplifier



a. Remove the curly cord from the phone base, and then connect it to the **HANDSET** port on the bottom of the hub.

b. Connect one end of the AC power adapter assembly to the **DC 9V** port on the left side of the hub, and then connect the other end of the assembly to an AC power outlet in the wall.



c. Connect the (short) Hub-to-base cable to the phone base port where the curly cord was connected. Connect the other end of the cable to the **BASE** port on the top of the hub.

If your hub is already on a compatible channel, the green light on the hub will be illuminated.

Corded Phone Setup

Check for Compatibility

- a. Pick up the telephone handset, and listen for a dial tone.
- b. Take one of these steps:

•



If you hear a dial tone: On the hub, slide the **VOLUME** adjuster until you hear a noticeable change in volume, or slide the **TONE** adjuster until you hear a noticeable change in tone quality. If volume or tone do not change, and the green light does not turn on, see below.

•



If you do not hear an adjustable dial tone: On the back of the hub, turn the compatibility dial very slowly until you hear a dial tone. On the front of the hub, slide the **VOLUME** adjuster until you hear a noticeable change in volume, or slide the **TONE** adjuster until you hear a noticeable change in tone quality.

Corded Phone Setup

Bluetooth Pairing and Disconnection

The first time that you use the QH2 system, you must pair the hub with the Bluetooth headset, streamer or neckloop that you will wear while using the system.

To pair the QH2 hub with a Bluetooth headset, streamer or neckloop

If necessary, consult your Bluetooth listening device user's guide to locate its multi-function button for pairing and to identify its pairing indicators. **Make sure there are no other Bluetooth devices nearby that have been paired to your headset, neckloop or streamer (Example: Cell Phone).** There should be no Bluetooth devices nearby that may interfere with your pairing and connection.

1. Wear your Bluetooth listening device, and remain within three feet of the QH2 hub.
2. On the hub, press and hold  until the blue and red light flashes.
3. Turn on your listening device, and then put it in the Bluetooth pairing mode. When the hub and listening device are paired, **the hub flashes a blue light every three seconds.** Your listening device also indicates when pairing is successful.
4. Confirm that you hear a dial tone by taking one of these steps **(Continued on next page):**

Corded Phone Setup

(For systems using a corded-phone with no headset mode and no installed handset lifter):

- a. On the hub, press and release .
- b. Manually lift the handset from the base. You should hear the dial tone in your Bluetooth listening device.
- c. Press and release . The dial tone stops.
- d. Manually return the handset to the base.

Note: There is no volume adjustment during a Bluetooth connection. Volume adjustments can be made with your Bluetooth headset, streamer, or neckloop.

(For all other corded system configurations):

- a. On the hub, press and release . You should hear the dial tone in your Bluetooth listening device.
- b. Press and release . The dial tone stops.

Corded Phone Setup - Answer a Call

Answer or place a call using a corded phone system

Note: Besides the following instructions for using the system with Bluetooth, you can simply pick up the telephone handset to answer and place calls. The hub only acts as an amplifier in this case, so that you can further amplify or adjust the volume and adjust the tone for the call by using the sliders on the hub.

To answer a call using a corded phone

1. Within 30 feet of the hub, put on your Bluetooth listening device.
2. On the hub, press and release .

Note: You must manually lift the handset from the phone base.

3. (Optional) Adjust the tone and volume:
 - Slide the TONE slider on the hub until you hear the tone that you prefer.
 - Adjust volume on your Bluetooth listening device or phone.

Note: If you increase volume on the listening device and the phone, and sound distortion occurs, lower the volume on the listening device or the phone.

4. To end the call, press and release  on the hub.

Note: You must manually return the handset to the base.

Corded Phone Setup - Place a Call

To place a call using a corded phone

1. Within 30 feet of the hub, put on your Bluetooth listening device.

2. On the hub, press and release .

You should hear the dial tone in your listening device.

Note: You must manually lift the handset from the phone base.

3. Press the desired keys on your telephone keypad.

4. (Optional) Adjust the tone and volume:

- Slide the TONE slider on the hub until you hear the tone that you prefer.
- Adjust volume on your Bluetooth listening device or phone.

Note: If you increase volume on the listening device and the phone, and sound distortion occurs, lower the volume on the listening device or the phone.

5. To end the call, press and release .

Note: You must manually return the handset to the base.

Corded Phone Setup - Break Bluetooth Connection

To break the Bluetooth connection between the hub and the listening device

If you find it necessary to break the Bluetooth connection, press and hold  on the hub until the red light flashes.

Otherwise, you should keep the Bluetooth power on, so that you are ready to place or answer calls.

You can answer or place a call when the Bluetooth listening device is within 30 feet of the hub.

To remain ready to answer and place calls, keep these components in this state:

- Bluetooth listening device (headset, streamer or neckloop) on and Bluetooth-connected
- QH2 hub connected to AC power and Bluetooth-connected

Note: Instead of pressing  on the hub to answer or end a call, you might be able to use a multi-function button on your Bluetooth listening device. For example, if you are using a ClearSounds Quattro neckloop, you can press the multifunction button twice to engage a call and once to disconnect a call. Check your listening device user's manual.

Cordless Phone Setup



Cordless Phone Setup - At Home

This setup is **only intended** for use with a single-line home **phone system**.

Connect the Hub

- a. Connect one end of the AC power adapter assembly to the DC 9V port on the left side hub, and then connect the other end of the assembly to an AC power outlet in the wall.
- b. Remove the telephones RJ11 jack from your wall phone line socket. Plug in the provided phone splitter.



- c. Connect the hub-to-splitter cable: Connect one end of the cable into the **LINE** port on the hub, and then connect the other end of the cable to a port on the phone line splitter.

The green light on the hub will not illuminate until you have connected a bluetooth headset, neckloop or hearing aid streamer to it and engaged the line. See page 17 for “Bluetooth Pairing and disconnection.”

- d. Pick up the telephone handset, and confirm that you hear a dial tone.

TIP: Make sure the compatibility dial on the back of the hub is turned to face the white dot for best results.

Cordless Phone Setup

Bluetooth pairing and disconnection

The first time that you use the QH2 system, you must pair the hub with the Bluetooth headset, streamer or neckloop that you will wear while using the system.

To pair the QH2 hub with a Bluetooth headset, streamer or neckloop

If necessary, consult your Bluetooth listening device user's guide to locate its multi-function button for pairing and to identify its pairing indicators.

1. Wear your Bluetooth listening device, and remain within three feet of the QH2 hub.
2. On the hub, press and hold  until the blue and red light flashes.
3. Turn on your listening device, and then put it in the Bluetooth pairing mode.

When the hub and listening device are paired, the hub flashes a blue light every three seconds. Your listening device also indicates when pairing is successful.

4. Confirm that you hear a dial tone by taking one of these steps:
 - a. On the hub, press and release . You should hear the dial tone in your Bluetooth listening device.
 - b. Manually lift the handset from the base.

Cordless Phone Setup - Answer a Call

To answer a call using a cordless phone

1. Within 30 feet of the hub, put on your Bluetooth listening device.
2. On the hub, press and release . You will hear the caller in your listening device.
3. **(Optional)** Adjust the tone and volume:
 - Slide the **TONE** slider on the hub until you hear the tone that you prefer.
4. Adjust volume on your Bluetooth listening device.

Note: If you increase volume on the listening device, and sound distortion occurs, lower the volume on the listening device the call, press and release .

5. To end the call, press and release .

Cordless Phone Setup - Place a Call

To place a call using a cordless phone

1. Within 30 feet of the hub, put on your Bluetooth listening device.
2. Lift your cordless handset, press the handset call button, and then press the desired keys on your telephone keypad.
3. When you hear the phone dialing through the handset, press and release  on the hub. When you hear the dialing in your listening device, the hub has transferred the call from the handset to your listening device.
4. Press the End button on your cordless phone, and return the handset to the base.
5. (Optional) Adjust the tone and volume:
 - Slide the **TONE** slider on the hub until you hear the tone that you prefer.
 - Adjust volume on your Bluetooth listening device.

Note: If you increase volume on the listening device, and sound distortion occurs, lower the volume on the listening device.
6. To end the call, press and release .

Cordless Phone Setup - Break Bluetooth Connection

To break the Bluetooth connection between the hub and the listening device

If you find it necessary to break the Bluetooth connection, press and hold  on the hub until the red light flashes.

Otherwise, you should keep the Bluetooth power on, so that you are ready to place or answer calls.

You can answer or place a call when the Bluetooth listening device is within 30 feet of the hub.

To remain ready to answer and place calls, keep these components in this state:

- Bluetooth listening device (headset, streamer or neckloop) on and Bluetooth-connected
- QH2 hub connected to AC power and Bluetooth-connected

Note: Instead of pressing  on the hub to answer or end a call, you might be able to use a multi-function button on your Bluetooth listening device. For example, if you are using a ClearSounds Quattro 4.0 neckloop, you can press the multifunction button twice to engage a call and once to disconnect a call. Check your listening device user's manual.

Troubleshooting

Symptom	Solution
The Bluetooth connection is disrupted (between the hub and the Bluetooth listening device).	Ensure the following: <ul style="list-style-type: none">• Your listening device is within 30 feet of the hub.• The listening device Bluetooth feature is on, and it is paired with the hub (the hub flashes a blue light every three seconds). See page 11.• The hub is connected to AC power and the hub's green light is on.
The sound is muffled.	If the volume is turned up on the phone and on the Bluetooth listening device, the sound can be muffled. Turn down the volume on the phone or listening device.
Sometimes the Bluetooth listening device is not charged, or it is not near the phone. Can I use the system without the Bluetooth listening device?	Yes. You can simply lift the phone handset and dial or speak. You can use the hub to adjust the volume and tone.

Compliance

Hearing Aid Compatibility (HAC)

This phone meets the U.S. Federal Communications Commission (FCC) standards for hearing aid compatibility.

FCC wants you to know

1. This equipment complies with Part 68 of the FCC rules and the requirements adopted by the ACTA. On the bottom of this equipment is a label that contains, among other information, a product identifier of [US: CLSTE12B]. If requested, this number must be provided to the phone company.
2. If the ClearSounds® QH2 Bluetooth Phone Hub System causes harm to the phone network, the phone company will notify you in advance that temporary discontinuance of service may be required. However, if advance notice is not practical, the phone company will notify the customer as soon as possible. Also, you will be advised of your right to file a complaint with the FCC if you believe it is necessary.
3. The phone company may make changes in its facilities, equipment, operations or procedures that could affect the operation of the equipment. If this happens, the phone company will provide advance notice in order for you to make necessary modification to maintain uninterrupted service.
4. If you experience trouble with this equipment, you should disconnect it from the network until the problem has been corrected or until you are sure that the equipment is not malfunctioning.
5. Please follow instructions for repairing if any (e.g. battery replacement section); otherwise do not alter or repair any parts of device except as specified.

Compliance

6. Connection to party line service is subject to state tariffs. Contact the state public utility commission, public service commission or corporation commission for information.
7. If the phone company requests information on what equipment is connected to their lines, inform them:
 - a. The phone number that this unit is connected to
 - b. The ringer equivalence number [1.2b]
 - c. The USOC jack required [RJ11C]
 - d. The FCC Registration Number [US:]

Items (b) and (d) are indicated on the label. The ringer equivalence number (REN) is used to determine how many devices can be connected to your phone line. In most areas, the sum of the RENs of all devices on any one line should not exceed five (5.0). If too many devices are attached, they may not ring properly.

Service Requirements

In the event of equipment malfunction, all repairs should be performed by our Company or an authorized agent. It is the responsibility of users requiring service to report the need for service to our Company or to one of our authorized agents. Service can be facilitated through our office at:

ClearSounds Communications, Inc.
1743 Quincy Avenue, Suite 143
Naperville, IL 60540

Canada Statement

Industry Canada CS03 Statement

This product meets the applicable Industry Canada technical specification. This is confirmed by the registration number. The abbreviation IC before the registration number signifies that registration was performed based on a Declaration of Conformity indicating that Industry Canada technical specifications were met. It does not imply that Industry Canada approved the equipment.

Before installing this equipment, users should ensure that it is permissible to be connected to the facilities of the local telecommunications company. The equipment must also be installed using an acceptable method of connections. The customer should be aware that compliance with the above conditions may not prevent degradation of service in some situations.

Repairs to certified equipment should be made by an authorized Canadian maintenance facility designated by the supplier. Any repairs or alterations made by the user of this equipment, or equipment malfunctions, may give the telecommunication company cause to request the user to disconnect the equipment.



WARNING: Users should not attempt to make such connections themselves, but should contact the appropriate electric inspection authority, or an electrician, as appropriate.



CAUTION: The Ringer Equivalence Number (REN) assigned to each terminal device provides an indication of the maximum number of terminals allowed to be connected to a phone interface.

The termination on an interface may consist of any combination of devices subject only to the requirement that the sum of the Ringer Equivalence Numbers of all the devices does not exceed 5. An alternative indication of the REN is the Load Number (LN); this number should be divided by 20 to approximate the equivalent REN. The REN Number of this unit is indicated on the bottom of your hub.

If your home has specially-wired alarm equipment connected to the phone line, ensure the installation of the ClearSounds QH2 Bluetooth Phone Hub System does not disable alarm equipment. If you have questions about what will disable alarm equipment, consult your phone provider or a qualified installer.

This Class B digital apparatus complies with Canadian ICES-003. The product is warranted by ClearSounds against manufacturing defects in material and workmanship under normal use for one (1) year from the date of purchase. Should you experience a problem, contact our Customer Service department or visit our website at www.clearsounds.com.

Specifications

Component	Specification
Ringer-detecting RDI port	15 Hz - 68 Hz
Bluetooth version	3.0
Frequency	2.4 G SIM, Channel 79
Modulation system	FSK
Power	Class 2, -6~4DBM Indoor conversion distance up to 15 meters, with only the PCB antenna (without the extra PA and extended antenna)

Warranty

Be sure to save your sales receipt as the proof of purchase date should you need warranty service. The serial number for this product is on the underside of the hub.

Your ClearSounds QH2 Bluetooth Phone Hub System comes with a (1) year limited warranty from the date of purchase. This products is warranted by ClearSounds against manufacturing defect in material and workmanship under normal use. In the event that this product fails to function properly due to defects in materials or workmanship within one year of the original purchase, return the unit (freight prepaid) with proof of purchase (sales receipt or packing slip, no exceptions) to ClearSounds Communications. ClearSounds Communications will either repair or replace the unit (with a refurbished unit or unit of equal condition) and return it you (using UPS/USPS ground shipping) at no cost to you. The warranty replacement or repair will be warranted for a period of 90 days or the remainder of the original warranty period (from original purchase date), whichever is longer.

The warranty does not cover accidents, negligence or breakage to any parts. This includes shipping damage, failure to follow instructions, misuse, fire, floods, use of incompatible accessories, acts of God or failure in your phone service carrier's line service. The product must not be tampered with or taken apart by anyone who is not an authorized ClearSounds Communications® representative. Tampering with the phone will void any written or implied warranties.

If a defect covered by this warranty should occur, promptly contact a Customer Service representative by phone, or log on to www.clearsounds.com to obtain a Return Merchandise Authorization (RMA) number and shipping instructions before shipping the product to us. Any shipment without an RMA number will not be accepted and will be returned to you at your expense. All authorized returned products must be accompanied with proof of purchase and a brief explanation of the problem.

For out-of-warranty repairs and service, please contact our Customer Service department for instructions. The repair shall be warranted for 90 days.

IMPORTANT: YOUR RECEIPT IS PART OF YOUR WARRANTY AND MUST BE RETAINED AND PRODUCED IN THE EVENT OF A WARRANTY CLAIM.

Accessories

Please visit the ClearSounds website (www.clearsounds.com) or call ClearSounds (800-965-9043) for additional accessories.



ClearSounds QT4L



Smart Sound Earbuds (RS062)

...asking more of ourselves, our lives, our minds and hearts. We can't just sit there, become ACTIVE. With the amazing events that unfolded in recent years, each and every one of us has been IM... a need to become part of something bigger than ourselves. ... cultures have all begun spiraling in a new direction... a direc... th no less an end result than changing the world as we know... ever. Technology brought citizens of the earth closer together. The organic change movement was helped by these devices. Wor... nnections were made. Stories were told from exotic locations. ... l new media intermingled and became something entirely new, a... communicate faster and more potent than ever before. Interco... rough technology, we've grown into bolder and wiser variations... ourselves. And we can no longer afford to live stagnant and bor... mmunication is about connection. It's about living life a little... plified than you once did. If you're not connecting with the... l the world around you, life quickly shifts away from some... ter and becomes something sweet. As the senior generation

© 2017 ClearSounds CS-QH2UG Ver 3