

## CLA7V2 Amplified Power Neckloop Accessory



### Your “Master Key” to connectivity!

Are you an iPhone, iTouch or soon-to-be iPad user? The ClearSounds CLA7v2 is the must-have accessory! If you wear hearing aids, the CLA7v2 maximizes everything these new mobile devices offer.

The ClearSounds® CLA7V2 Amplified Neckloop is a convenient little device you wear around your neck and plug into a cell phone, cordless phone, MP3 player, laptop, gaming device or any audio device.

The handy accessory gives you adjustable amplification up to 30 dB and reduces or eliminates interference between cellular phones and hearing aid processors. It also features a built in microphone for hands-free conversation.

You need not purchase or switch between multiple headsets for different activities. Four connection cables are included to accommodate the most common types of phones and audio equipment. The amplified neckloop is powered by AAA batteries to provide approximately 200 hours of operation.



### KEY FEATURES

Mobile, home & office applications

30 dB adjustable amplification

User-friendly controls & cables

One device, multiple applications: cellular phones, home/office phones, laptop computers, portable gaming devices, iPods, personal/TV listening equipment and other audio devices

### Additional Features

- No need to purchase multiple headsets or switch audio headsets between different activities
- Approximately 200 hours of operation on 2 AAA batteries
- Reduces or eliminates interference between cellular phones and hearing aid(s) or cochlear processor(s)
- Requires T-coil(s)
- Best sound quality
- Includes 4 connection cables
- T-coil compatible
- Binaural design
- Push-button slider to cinch neckloop for closer microphone placement
- Trouble-free operation
- Offers a 30 dB max gain
- Rotary volume control for easy adjustment and multiple levels of amplification
- Built-in microphone
- Quick release feature
- Battery powered with (2) AAA batteries (included)
- Size: 1.25" x 2.5" x .75"
- Color: Black
- 1-year warranty

### Great Value! Four Connection Cables

- Each cable is removeable from the base of the CLA7V2
- 3.5mm to 2.5mm mono hands-free cable – used for hands-free connections with cordless and cellular telephones with a mono 2.5mm jack (phone may have an optional 2.5mm mono adapter if the jack is not built-in).
- 3.5mm to 2.5mm stereo hands-free cable – used with newer cellular phones that have stereo 2.5mm jacks (phone may have an optional 2.5mm stereo adapter if the jack is not built-in).
- 3.5mm to 3.5mm mono audio cable – used with corded phones that provide a 3.5mm mono jack for connection to neckloops and headsets.
- 3.5mm to 3.5mm stereo audio/hands-free cable – used for streaming music and voice, directly to your iPhone, Blackberry or other audio device.