

Total Wireless Lifestyle Connectivity

The Quattro Amplified Bluetooth Neckloop and QLink Accessories provide an unmatched system for ClearSounds Total 360 LifeStyle Connectivity™ that is now useable by everyone



CLEAROUNDS.COM

PC

Listen to audio and make VoIP phone calls using SKYPE

Parts Checklist:

- QUATTRO Amplified Bluetooth Neckloop
 - RS062 SmartSound Earbuds (for non hearing aid, with t-coil, users)
 - QLink Stereo Bluetooth Transmitter (please note that you will need to use a webcam with built in microphone)
 - Bluetooth USB Micro Adapter (Optional)
- To use the microphone in the QUATTRO and not in the Webcam

LAPTOP

Listen to audio and make VoIP phone calls using SKYPE

Parts Checklist:

- QUATTRO Amplified Bluetooth Neckloop
 - RS062 SmartSound Earbuds (for non hearing aid, with t-coil, users)
 - QLink Stereo Bluetooth Transmitter (please note that you will need to use a webcam with built in microphone)
 - Bluetooth USB Micro Adapter (Optional)
- To use the microphone in the QUATTRO and not in the Webcam



CELLULAR

Wireless handsfree conversation. Listen to the music on your Smart phone. Works with all Bluetooth Enabled Cellular phones.

Parts Checklist:

- QUATTRO Amplified Bluetooth Neckloop
- RS062 SmartSound Earbuds (for non hearing aid, with t-coil, users)

TV LISTENING

Works with all TV's including LCD, LED and Plasma. (TV and cell phone can be connected at the same time)

Parts Checklist:

- QUATTRO Amplified Bluetooth Neckloop
- RS062 SmartSound Earbuds (for non hearing aid, with t-coil, users)
- QLink Stereo Bluetooth Transmitter

iPOD/AUDIO/MP3

Smart technology allows you to connect your MP3 player and cellular handset at the same time. Your music will mute to answer an incoming call.

Parts Checklist:

- QUATTRO Amplified Bluetooth Neckloop
- RS062 SmartSound Earbuds (for non hearing aid, with t-coil, users)
- QLink Stereo Bluetooth Transmitter (only required if your MP3 player does not have Bluetooth)