

QH2 BLUETOOTH HUB & PHONE AMPLIFIER

With the QH2, you can connect a Bluetooth headset, neckloop or hearing aid streamer to your home or office phone for hands-free and cord-free conversation. The QH2 is compatible with virtually all corded desk phones, as well as corded and cordless home phones. Various set-up options provide superior flexibility.

The QH2 provides 30 dB of adjustable amplification and frequency management in both the phone handset and Bluetooth headset so you can be comfortable knowing amplification is always available. Easily switch calls between the phone handset and Bluetooth headset.

* Bluetooth headset, neckloop or streamer must support Hands-Free Profile



FEATURES

- Bluetooth 3.0
- AC powered
- 33ft wireless range
- Handset lifter included
- Stores up to 8 paired devices
- Remarkable sound quality

- Supports remote answering**
- Pairs with any Bluetooth headset
- 30dB amplification & frequency control
- Stylish design with a 'soft touch' finish

**Depending on your Bluetooth headset

Clear

QH2 BLUETOOTH HUB & PHONE AMPLIFIER

The hub will convert any corded or cordless residential telephone into a Bluetooth phone, so you can pair it with a compatible Bluetooth Headset, streamer or neckloop for wireless communication in your home.



The hub will also convert any corded business (PBX, VoIP etc.) telephone into a Bluetooth phone so you can pair it with a compatible Bluetooth Headset for wireless communication in your office.

Get traditional amplification for your office or home corded telephone system while enjoying the benefits of Bluetooth wireless connectivity.